

---

## Product Datasheet

### Human IL-17A ELISpot Kit 856.151PC

<b>Product Name</b>	Human IL-17A ELISpot Kit
<b>Catalog Number</b>	856.151PC
<b>Description</b>	Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Specificity</b>	Recognizes natural human IL-17A