

Product Datasheet

Human IL-17A/F ELISpot Set 856.171.005S

Product Name	Human IL-17A/F ELISpot Set
Catalog Number	856.171.005S
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 5 sterile plates.
Tested Applications	ELISpot
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Note	Spots read visually or with reader
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes natural human IL-17A/F