

Product Datasheet

Anti-h E2 (Estradiol) 10450 700011

Product Name	Anti-h E2 (Estradiol) 10450
Catalog Number	700011
Description	Monoclonal antibody raised against human 17b Estradiol. In vitro cultured mouse monoclonal antibody.
Species/Host	Mouse
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	> 1.0 mg/ml
Storage	+2-8 °C
Isotype	IgG1
Clonality	Monoclonal
Epitope	Not determined
Affinity constant	N/D
Buffer	0.1 M PBS, pH 7.4, 0.9 % NaCl, 0.05 % NaN3 as a preservative
Cross Reactivity	Not determined
Specificity	Antibody recognizes human estradiol (E2)
Shelf Life	10 years