

## Product Datasheet

### Anti-h Insulin 9351 700014

<b>Product Name</b>	Anti-h Insulin 9351
<b>Catalog Number</b>	700014
<b>Description</b>	Monoclonal antibody raised against human Insulin. In vitro cultured mouse monoclonal antibody.
<b>Tested Applications</b>	ELISA
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	> 1.0 mg/ml
<b>Storage</b>	+2-8°C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Epitope</b>	N/D
<b>Purity</b>	≥ 95 %
<b>Buffer</b>	0.1 M PBS, pH 7.4, 0.9 % NaCl, 0.05 % NaN <sub>3</sub> as a preservative
<b>IEF Profile</b>	6.8 - 7.2
<b>Cross Reactivity</b>	N/D
<b>Specificity</b>	Antibody recognizes human insulin
<b>Matched Pairs</b>	Anti-h Insulin 9350 (700013)
<b>Shelf Life</b>	10 years