

# Product Datasheet

## Rat TNF alpha ELISpot Set 867.021.001P

Product Name	Rat TNF alpha ELISpot Set
Catalog Number	867.021.001P
Description	Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 1 non-sterile plate.
Tested Applications	ELISpot
Alternative Names	TNFa
Brand	Diaclone
Storage	Store at +4°C
Note	Spots read visually or with reader
Cross Reactivity	No cross reactivity with other rat cytokines
Specificity	Recognizes natural rat TNF-a