

Product Datasheet

Rat TNF alpha ELISpot Kit 867.021PC

Product Name	Rat TNF alpha ELISpot Kit
Catalog Number	867.021PC
Description	Diaclone Pre-coated ELISpot kits include precoated PVDF plates, Detection antibody, Alkaline phosphatase conjugate, BSA, BCIP/NBT ready-to-use substrate buffer.
Tested Applications	ELISpot
Alternative Names	TNFa
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other rat cytokines
Specificity	Recognizes natural rat TNF-a