

Product Datasheet

Rat TNF alpha ELISpot Set 867.021.001S

| | |
|----------------------------|---|
| Product Name | Rat TNF alpha ELISpot Set |
| Catalog Number | 867.021.001S |
| Description | Diaclone ELISpot Sets include capture and detection antibodies, Streptavidin - Alkaline Phosphatase conjugated, BSA, BCIP/NTB, blocking reagent, and 1 sterile plate. |
| Tested Applications | ELISpot |
| Alternative Names | TNFa |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Note | Spots read visually or with reader |
| Cross Reactivity | No cross reactivity with other rat cytokines |
| Specificity | Recognizes natural rat TNF-a |