

## Product Datasheet

### Human IFN $\gamma$ / IL-10 Dual ELISpot set 874.000.015S

<b>Product Name</b>	Human IFN $\gamma$ / IL-10 Dual ELISpot set
<b>Catalog Number</b>	874.000.015S
<b>Description</b>	Diaclone Dual ELISpot Sets include capture and detection antibodies for both of the targeted cytokines, Streptavidin-Alkaline Phosphatase conjugated, anti-FITC HRP conjugated mAb, BSA, BCIP/NTB and AEC, blocking reagent and 15 sterile plates.
<b>Tested Applications</b>	ELISpot
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity with other human cytokines. Cross reactivity with simian IFN-g. Cross reactivity with simian IL-10. No cross reactivity with viral IL-10
<b>Specificity</b>	Recognizes natural human IL-10 and IFN-g