

Product Datasheet

Human IFN γ / IL-4 Dual ELISpot Set 874.010.005P

Product Name	Human IFN γ / IL-4 Dual ELISpot Set
Catalog Number	874.010.005P
Description	Diaclone Dual ELISpot Sets include capture and detection antibodies for both of the targeted cytokines, Streptavidin-Alkaline Phosphatase conjugated, anti-FITC HRP conjugated mAb, BSA, BCIP/NTB and AEC, Blocking reagent and 5 non-sterile plates.
Tested Applications	ELISpot
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IFN-g
Specificity	Recognizes natural human IL-4 and IFN-g