

Product Datasheet

Anti-h PTH 10750 700024

Product Name	Anti-h PTH 10750
Catalog Number	700024
Description	Monoclonal antibody raised against human Parathyroid Hormone (PTH). In vitro cultured mouse monoclonal antibody.
Species/Host	Mouse
Tested Applications	RIA
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	> 1.0 mg/ml
Storage	+2-8 °C
Isotype	IgG1
Clonality	Monoclonal
Epitope	Amino acids 57-61
Affinity constant	N/D
Buffer	0.1 M PBS, pH 7.4, 0.9 % NaCl, 0.05 % NaN ₃ as a preservative
IEF Profile	N/D
Cross Reactivity	Not determined
Specificity	Antibody recognizes human parathyroid hormone (PTH)
Shelf Life	10 years