

Product Datasheet

Human IFN γ / IL-5 Dual ELISpot Set 874.020.020S

Product Name	Human IFN γ / IL-5 Dual ELISpot Set
Catalog Number	874.020.020S
Description	Diaclone Dual ELISpot Sets include capture and detection antibodies for both of the targeted cytokines, Streptavidin-Alkaline Phosphatase conjugated, anti-FITC HRP conjugated mAb, BSA, BCIP/NTB and AEC, blocking reagent and 20 sterile plates.
Tested Applications	ELISpot
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	Cross-reactive with simian IFN- γ
Specificity	Recognizes both natural and recombinant human IFN- γ and IL-5