

Product Datasheet

Human IL-10 / IL-4 Dual Fluorospot Set 874.022.015P

Product Name	Human IL-10 / IL-4 Dual Fluorospot Set
Catalog Number	874.022.015P
Description	Diaclone Dual Fluorospot Sets include capture antibody for cytokine 1, FITC-conjugated detection antibody for cytokine 1, anti-FITC antibody green fluorescence conjugate, capture antibody for cytokine 2, biotinylated detection antibody for cytokine 2, streptavidin-phycoerythrin conjugate, BSA. Kit can be sold with or without plates.
Tested Applications	Fluorospot
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant human IL-4 and IL-10