

Product Datasheet

Human IL-10 / IL-4 Dual ELISpot Set 874.030.001P

Product Name	Human IL-10 / IL-4 Dual ELISpot Set
Catalog Number	874.030.001P
Description	Diaclone Dual ELISpot Sets include capture and detection antibodies for both of the targeted cytokines, Streptavidin-Alkaline Phosphatase conjugated, anti-FITC HRP conjugated mAb, BSA, BCIP/NTB and AEC, Blocking reagent and 1 non-sterile plate.
Tested Applications	ELISpot
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant human IL-4 and IL-10