

Product Datasheet

Human IFN γ / IL-5 Dual Fluorospot Set 874.032.001P

Product Name	Human IFN γ / IL-5 Dual Fluorospot Set
Catalog Number	874.032.001P
Description	Diaclone Dual Fluorospot Sets include capture antibody for cytokine 1, FITC-conjugated detection antibody for cytokine 1, anti-FITC antibody green fluorescence conjugate, capture antibody for cytokine 2, biotinylated detection antibody for cytokine 2, streptavidin-phycoerythrin conjugate, BSA. Kit can be sold with or without plates.
Tested Applications	Fluorospot
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	Cross-reactive with simian IFN-g
Specificity	Recognizes both natural and recombinant human IFN-g and IL-5