

Product Datasheet

Anti-h CEA 5912 SP-1 100049

Product Name	Anti-h CEA 5912 SP-1
Catalog Number	100049
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	CLIA, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	1.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Note	In Bhayana et al. (1989) authors made an epitope mapping for 11 monoclonal antibodies from Medix Biochemica and they found that the antibodies bound to 5 different epitope groups (A to F) and 5912 was found to bind to epitope group F.
Isotype	IgG1
Clonality	Monoclonal
Epitope	Epitope group F as described in Bhayana et al. (1989)
Purity	≥ 95 %
Affinity constant	KA= 4 x 10 ⁹ 1/M
Associated Products	Native CEA antigen, Lee Biosolutions 151-09 and 151-10
Buffer	0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.8-7.6
Cross Reactivity	Recognizes human NCA-2 (non-specific cross-reacting antigen-2)
Specificity	Antibody recognizes human carcinoembryonic antigen
Shelf Life	36 months
References	Bhayana, V. and Diamandis, E.P. (1989) A double monoclonal time-resolved immunofluorometric assay of carcinoembryonic antigen in serum. Clin. Biochem., 22:433-438