

## Product Datasheet

### Human IFN $\gamma$ / IL-2 Dual ELISpot Set 874.040.001P

<b>Product Name</b>	Human IFN $\gamma$ / IL-2 Dual ELISpot Set
<b>Catalog Number</b>	874.040.001P
<b>Description</b>	Diaclone Dual ELISpot Sets include capture and detection antibodies for both of the targeted cytokines, Streptavidin-Alkaline Phosphatase conjugated, anti-FITC HRP conjugated mAb, BSA, BCIP/NTB and AEC, Blocking reagent and 1 non-sterile plate.
<b>Tested Applications</b>	ELISpot
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	Cross reactivity with simian IFN-g
<b>Specificity</b>	Recognizes both natural and recombinant human IFN-g and IL-2