

# Product Datasheet

## Human IFN $\gamma$ / IL-2 Dual ELISpot Set 874.040.005P

Product Name	Human IFN $\gamma$ / IL-2 Dual ELISpot Set
Catalog Number	874.040.005P
Description	Diaclone Dual ELISpot Sets include capture and detection antibodies for both of the targeted cytokines, Streptavidin-Alkaline Phosphatase conjugated, anti-FITC HRP conjugated mAb, BSA, BCIP/NTB and AEC, Blocking reagent and 5 non-sterile plates.
Tested Applications	ELISpot
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	Cross reactivity with simian IFN-g
Specificity	Recognizes both natural and recombinant human IFN-g and IL-2