

---

## Product Datasheet

### Human IL-10 / IL-2 Dual ELISpot Set 874.050.010P

<b>Product Name</b>	Human IL-10 / IL-2 Dual ELISpot Set
<b>Catalog Number</b>	874.050.010P
<b>Description</b>	Diaclone Dual ELISpot Sets include capture and detection antibodies for both of the targeted cytokines, Streptavidin-Alkaline Phosphatase conjugated, anti-FITC HRP conjugated mAb, BSA, BCIP/NTB and AEC, blocking reagent and 10 non-sterile plates.
<b>Tested Applications</b>	ELISpot
<b>Alternative Names</b>	NULL
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Specificity</b>	Recognizes both natural and recombinant human IL-10 and IL-2