

## Product Datasheet

### Recombinant Lactoferrin, 50 µg 710016

<b>Product Name</b>	Recombinant Lactoferrin, 50 µg
<b>Catalog Number</b>	710016
<b>Description</b>	Recombinant human lactoferrin, also known as lactotransferrin/LTF antigen is produced by mammalian expression system and the target gene encoding Gly20-Lys710 is expressed with a 6His tag at the C-terminus.
<b>Species/Host</b>	Human
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Lyophilized
<b>Concentration</b>	0.5 mg/ml when reconstituted with 100 µl of deionized water
<b>Storage</b>	-20 °C
<b>Purity</b>	SDS-PAGE > 90 %
<b>Buffer</b>	Lyophilized from 20 mM Na-phosphate, pH 7.4, 150 mM NaCl
<b>Molecular Weight</b>	Predicted molecular weight: 77 kDa
<b>Shelf Life</b>	Unspecified