

Product Datasheet

Human IL-17A/IL-4 Dual Fluorospot Set 874.092.001P

Product Name	Human IL-17A/IL-4 Dual Fluorospot Set
Catalog Number	874.092.001P
Description	Diaclone Dual Fluorospot Sets include capture antibody for cytokine 1, FITC-conjugated detection antibody for cytokine 1, anti-FITC antibody green fluorescence conjugate, capture antibody for cytokine 2, biotinylated detection antibody for cytokine 2, streptavidin-phycoerythrin conjugate, BSA. Kit can be sold with or without plates.
Tested Applications	Fluorospot
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines.
Specificity	Recognizes natural human IL-17A and IL-4