

Product Datasheet

Recombinant Adiponectin, 500 µg 710031

Product Name	Recombinant Adiponectin, 500 µg
Catalog Number	710031
Description	Recombinant human adiponectin antigen is produced by mammalian expression system and the target gene encoding Glu19-Asn244 is expressed with a 6His tag at the C-terminus.
Species/Host	Human
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Concentration	0.5 mg/ml when reconstituted with 1000 µl of deionized water
Storage	-20 °C
Purity	SDS-PAGE > 90 %
Associated Products	Anti-h Adiponectin 1901 Anti-h Adiponectin 1902 Anti-h Adiponectin 1903
Buffer	Lyophilized from 20 mM Na-phosphate, pH 7.2, 150 mM NaCl
Molecular Weight	Predicted molecular weight: 30 kDa
Shelf Life	Unspecified