

Product Datasheet

Recombinant ApoA1, 1000 µg 710032

Product Name	Recombinant ApoA1, 1000 µg
Catalog Number	710032
Description	Recombinant human apolipoprotein A1 antigen is produced by mammalian expression system and the target gene encoding Arg19-Gln267 is expressed with a 6His tag at the C-terminus.
Species/Host	Human
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Concentration	0.5 mg/ml when reconstituted with 2000 µl of deionized water
Storage	-20 °C
Purity	SDS-PAGE > 90 %
Associated Products	Anti-h ApoA1 2001 Anti-h ApoA1 2002
Buffer	Lyophilized from 20 mM Na-phosphate, pH 7.2, 150 mM NaCl
Molecular Weight	Predicted molecular weight: 35 kDa
Shelf Life	Unspecified