

Product Datasheet

Recombinant Calprotectin S100A8, 500 µg 710035

Product Name	Recombinant Calprotectin S100A8, 500 µg
Catalog Number	710035
Description	Recombinant human S100 calcium binding protein A8 antigen is produced by E.coli expression system and the target gene encoding Met1-Glu93 is expressed with a 6His tag at the C-terminus.
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Concentration	0.5 mg/ml when reconstituted with 1000 µl of deionized water
Storage	-20 °C
Purity	SDS-PAGE > 90 %
Associated Products	Anti-h Calprotectin 3404 Anti-h Calprotectin 3405 Anti-h Calprotectin 3406
Buffer	Lyophilized from 20 mM HEPES, pH 7.4, 150 mM NaCl
Molecular Weight	Predicted molecular weight: 14 kDa
Shelf Life	Unspecified