

Product Datasheet

Recombinant Calprotectin S100A9, 1000 µg 710036

Product Name	Recombinant Calprotectin S100A9, 1000 µg
Catalog Number	710036
Description	Recombinant human S100 calcium binding protein A9 antigen is produced by E.coli expression system and the target gene encoding Thr2-Pro114 is expressed.
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Concentration	0.5 mg/ml when reconstituted with 2000 µl of deionized water
Storage	-20 °C
Purity	SDS-PAGE > 90 %
Associated Products	Anti-h Calprotectin 3403
Buffer	Lyophilized from PBS, pH 7.4
Molecular Weight	Predicted molecular weight: 13 kDa
Shelf Life	Unspecified