

Product Datasheet

Recombinant IL-6, 1000 µg 710038

Product Name	Recombinant IL-6, 1000 µg
Catalog Number	710038
Description	Recombinant human interleukin-6 antigen is produced by E.coli expression system and the target gene encoding Pro29-Met212 is expressed.
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Concentration	0.5 mg/ml when reconstituted with 2000 µl of deionized water
Storage	-20 °C
Purity	SDS-PAGE > 90 %
Associated Products	Anti-h IL-6 2703, Anti-h IL-6 2704, Anti-h IL-6 2706, Anti-h IL-6 2707
Buffer	Lyophilized from 20 mM Hac-NaAc, pH 5.0, 150 mM NaCl
Molecular Weight	Predicted molecular weight: 19 kDa
Shelf Life	Unspecified