

Product Datasheet

Recombinant Lactoferrin, 1000 µg 710040

Product Name	Recombinant Lactoferrin, 1000 µg
Catalog Number	710040
Description	Recombinant human lactoferrin, also known as lactotransferrin/LTF, antigen is produced by mammalian expression system and the target gene encoding Gly20-Lys710 is expressed with a 6His tag at the C-terminus.
Species/Host	Human
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Concentration	0.5 mg/ml when reconstituted with 2000 µl of deionized water
Storage	-20 °C
Purity	SDS-PAGE > 90 %
Buffer	Lyophilized from 20 mM Na-phosphate, pH 7.4, 150 mM NaCl
Molecular Weight	Predicted molecular weight: 77 kDa
Shelf Life	Unspecified