

Product Datasheet

Recombinant proBNP, 500 µg
710043

Product Name	Recombinant proBNP, 500 µg
Catalog Number	710043
Description	Recombinant human natriuretic peptide B antigen is produced by E. coli expression system and the target gene encoding His27-His134 is expressed with a 6His tag at the N-terminus. Predicted molecular weight: 15 kDa.
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Liquid
Concentration	0.5 - 1.5 mg/ml
Storage	-20 °C
Purity	SDS-PAGE > 90 %
Associated Products	Anti-h NT proBNP 1306 Anti-h NT proBNP 1307 Anti-h NT proBNP 1308 Anti-h NT proBNP 1309 Anti-h NT proBNP 1310 Anti-h NT proBNP 1311 Anti-h NT proBNP 1312
Buffer	20 mM Na-phosphate, pH 7.4, 150 mM NaCl, 1 mM EDTA, 1 mM DTT, 20 % Glycerol
Molecular Weight	Predicted molecular weight: 15 kDa
Shelf Life	Unspecified