

Product Datasheet

Recombinant S100B, 1000 µg 710044

Product Name	Recombinant S100B, 1000 µg
Catalog Number	710044
Description	Recombinant human S100 calcium binding protein B antigen is produced by E. coli expression system and the target gene encoding Met1-Glu92 is expressed with a 6His tag at the N-terminus.
Species/Host	E. coli
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Concentration	0.5 mg/ml when reconstituted with 2000 µl of deionized water
Storage	-20 °C
Purity	SDS-PAGE > 90 %
Associated Products	Anti-h S100B 11401 Anti-h S100B 11402 Anti-h S100B 11404
Buffer	Lyophilized from 20 mM PBS, pH 7.4, 150 mM NaCl
Molecular Weight	Predicted molecular weight: 12 kDa
Shelf Life	Unspecified