

# Product Datasheet

## Recombinant ST2, 1000 µg 710046

Product Name	Recombinant ST2, 1000 µg
Catalog Number	710046
Description	Recombinant human ST2, also known as interleukin-1 receptor-like 1, antigen is produced by mammalian expression system and the target gene encoding Lys19-Phe328 is expressed with a 6His tag at the C-terminus.
Species/Host	Human
Brand	Medix Biochemica
Form/Appearance	Lyophilized
Concentration	0.5 mg/ml when reconstituted with 2000 µl of deionized water
Storage	-20 °C
Purity	SDS-PAGE > 90 %
Associated Products	Anti-h ST2 10201 Anti-h ST2 10202 Anti-h ST2 10203 Anti-h ST2 10204 Anti-h ST2 10205 Anti-h ST2 10206 Anti-h ST2 10207
Buffer	Lyophilized from 20 mM Na-phosphate, pH 7.4, 150 mM NaCl
Molecular Weight	Predicted molecular weight: 60 kDa
Shelf Life	Unspecified