

## Product Datasheet

### Vitamin D antigen, Biotin conjug, 1 mg 710049

<b>Product Name</b>	Vitamin D antigen, Biotin conjug, 1 mg
<b>Catalog Number</b>	710049
<b>Description</b>	Vitamin D antigen, Biotin conjug, 1 mg
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Lyophilized antigen, on the bottom of the vial. Not visually detectable.
<b>Concentration</b>	1 mg/ml when reconstituted with 1000 µl of EtOH
<b>Storage</b>	Storage at -20 °C. After reconstitution 6 months at -20 °C. Long term storage at -80 °C.
<b>Source</b>	Synthetic
<b>Associated Products</b>	Anti-Vitamin 25 OH D2/D3 3850 Anti-Vitamin 25 OH D2/D3 3851
<b>Shelf Life</b>	60 months