

## Product Datasheet

### Recombinant SARS-CoV-2 N, 100 µg 710081

<b>Product Name</b>	Recombinant SARS-CoV-2 N, 100 µg
<b>Catalog Number</b>	710081
<b>Description</b>	Recombinant human coronavirus SARS-CoV-2 nucleocapsid protein with a C-terminal histidine tag.
<b>Species/Host</b>	CHO Cells
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.5-1 mg/ml
<b>Storage</b>	-70 °C
<b>Purity</b>	SDS-PAGE > 90 %
<b>Associated Products</b>	Anti SARS-CoV-2 N 3512 Anti SARS-CoV-2 N 3511 Anti SARS-CoV-2 N 3510 Anti SARS-CoV-2 N 3508 Anti SARS-CoV-2 N 3507
<b>Buffer</b>	PBS, pH 7.4
<b>Molecular Weight</b>	Predicted molecular weight: 47 kDa
<b>Shelf Life</b>	Unspecified