

Product Datasheet

Anti-Mouse IL-23 Monoclonal Antibody, Azide Free Clone B-B23 400.000

| | |
|----------------------------|--|
| Product Name | Anti-Mouse IL-23 Monoclonal Antibody, Azide Free Clone B-B23 |
| Catalog Number | 400.000 |
| Description | Mouse monoclonal |
| Species/Host | Mouse |
| Conjugation | Azide free |
| Tested Applications | ELISA |
| Immunogen | Recombinant murine IL-23p19 |
| Brand | Diaclone |
| Form/Appearance | Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free |
| Storage | Stable at +2-8°C for 12 months. For longer storage freeze aliquots. |
| Isotype | Mouse IgG1 Kappa light chain |
| Clonality | Monoclonal |
| Specificity | Recognises both natural and recombinant protein |
| Hybridoma | Myeloma X63/AG.8653 x C57 Black/6 spleen cells |