

## Product Datasheet

### Anti-Mouse IL-23 Monoclonal Antibody, Azide Free Clone B-B23 400.000

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Anti-Mouse IL-23 Monoclonal Antibody, Azide Free Clone B-B23                               |
| <b>Catalog Number</b>      | 400.000  |
| <b>Description</b>         | Mouse monoclonal   |
| <b>Species/Host</b>        | Mouse  |
| <b>Conjugation</b>         | Azide free   |
| <b>Tested Applications</b> | ELISA  |
| <b>Immunogen</b>           | Recombinant murine IL-23p19  |
| <b>Brand</b>               | Diaclone   |
| <b>Form/Appearance</b>     | Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free |
| <b>Storage</b>             | Stable at +2-8°C for 12 months. For longer storage freeze aliquots.                        |
| <b>Isotype</b>             | Mouse IgG1 Kappa light chain   |
| <b>Clonality</b>           | Monoclonal   |
| <b>Specificity</b>         | Recognises both natural and recombinant protein  |
| <b>Hybridoma</b>           | Myeloma X63/AG.8653 x C57 Black/6 spleen cells   |