

# Product Datasheet

## Human TGF beta 1 ELISA Kit 650.010

Product Name	Human TGF beta 1 ELISA Kit
Catalog Number	650.010
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.3-2000 pg/ml
Sensitivity	9 pg/ml
Tested Applications	ELISA
Alternative Names	TGFb1, TGF-beta1
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	Cross reactivity with sheep and pig TGF-β1
Specificity	Recognizes both natural and recombinant human TGF-β1
Sample Types	Serum, Plasma, Cell culture supernatant