

Product Datasheet

Human IL-22 ELISA Kit 650.140

Product Name	Human IL-22 ELISA Kit
Catalog Number	650.140
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.3-2000pg/ml
Sensitivity	5.0pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant human IL-22
Sample Types	Cell culture supernatants. Human serum and plasma.