

## Product Datasheet

### Human IL-22 ELISA Kit 650.140.096

<b>Product Name</b>	Human IL-22 ELISA Kit
<b>Catalog Number</b>	650.140.096
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	31.3-2000pg/ml
<b>Sensitivity</b>	5.0pg/ml
<b>Tested Applications</b>	ELISA
<b>Alternative Names</b>	NULL
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Specificity</b>	Recognizes both natural and recombinant human IL-22
<b>Sample Types</b>	Cell culture supernatants. Human serum and plasma.