

# Product Datasheet

## Murine IL-4 ELISA Kit 660.020

Product Name	Murine IL-4 ELISA Kit
Catalog Number	660.020
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	3.9-250 pg/ml
Sensitivity	3 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant murine IL-4
Sample Types	Serum, Plasma, Cell culture supernatant