

Product Datasheet

Murine Eotaxin ELISA Kit 660.060

Product Name	Murine Eotaxin ELISA Kit
Catalog Number	660.060
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-2000 pg/ml
Sensitivity	3.6pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant murine Eotaxin
Sample Types	Cell culture supernatant and serum tested, other body fluids suitable