

## Product Datasheet

### Murine Eotaxin ELISA Kit 660.060.096

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Murine Eotaxin ELISA Kit   |
| <b>Catalog Number</b>      | 660.060.096  |
| <b>Description</b>         | Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure. |
| <b>Range</b>               | 31.25-2000 pg/ml   |
| <b>Sensitivity</b>         | 3.6pg/ml   |
| <b>Tested Applications</b> | ELISA  |
| <b>Alternative Names</b>   | NULL   |
| <b>Brand</b>               | Diaclone   |
| <b>Storage</b>             | Store at +4°C  |
| <b>Specificity</b>         | Recognizes both natural and recombinant murine Eotaxin   |
| <b>Sample Types</b>        | Cell culture supernatant and serum tested, other body fluids suitable  |