Medix Biochemica

Oy Medix Biochemica Ab

Headquarters Klovinpellontie 3 FI-02180 Espoo

Finland

Email: sales@medixbiochemica.com

Medix Biochemica USA Inc.

10850 Metro Court Maryland Heights, MO 63043 United States of America

Email: sales@medixbiochemica.com

Product Datasheet

Murine CD40L ELISA Kit 660.070

Product Name Murine CD40L ELISA Kit

Catalog Number 660.070

Description

Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed

procedure.

Range

0.63-20ng/ml

Sensitivity

0.14ng/ml

Alternative

TRAP / CD154

Names

Diaclone

Storage

Brand

Store at +4°C

Specificity

Recognizes both natural and recombinant murine CD40L

Sample Types Cell culture supernatant and serum tested, other body fluids

suitable