

# Product Datasheet

## Murine CD40L ELISA Kit 660.070

Product Name	Murine CD40L ELISA Kit
Catalog Number	660.070
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	0.31-20ng/ml
Sensitivity	0.14ng/ml
Tested Applications	ELISA
Alternative Names	TRAP /CD154
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant murine CD40L
Sample Types	Cell culture supernatant and serum tested, other body fluids suitable