

Product Datasheet

Murine IL-13 ELISA Kit 660.080

| | |
|---------------------|--|
| Product Name | Murine IL-13 ELISA Kit |
| Catalog Number | 660.080 |
| Description | Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure. |
| Range | 2.8-500pg/ml |
| Sensitivity | 2.8pg/ml |
| Tested Applications | ELISA |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Specificity | Recognizes both natural and recombinant murine IL-13 |
| Sample Types | Cell culture supernatant and serum tested, other body fluids suitable |