

Product Datasheet

Murine IL-13 ELISA Kit 660.080.096

Product Name	Murine IL-13 ELISA Kit
Catalog Number	660.080.096
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	2.8-500pg/ml
Sensitivity	2.8pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant murine IL-13
Sample Types	Cell culture supernatant and serum tested, other body fluids suitable