

Product Datasheet

Goat Immunoglobulin G (IgG) H + L Antibody A151

Product Name Goat Immunoglobulin G (IgG) H + L Antibody

Catalog Number A151

Species/Host Rabbit

Tested Applications ELISA, WB, IEP

Immunogen Purified Goat IgG

Brand EastCoast Bio

Form/Appearance Affinity Purified Polyclonal Antibody

Storage 2-8°C. Do not freeze.

Clonality Polyclonal

Purity >95%

Buffer 10 mM Phosphate Buffered Saline, pH 7.2

Preservatives 0.1% Sodium Azide