

Product Datasheet

Vero Cells Antibody, FITC AG232

Product Name	Vero Cells Antibody, FITC
Catalog Number	AG232
Species/Host	Rabbit
Tested Applications	FC, IP, WB, IFA
Immunogen	Vero-76 cells (ATCC CRL-1587)
Brand	EastCoast Bio
Form/Appearance	Polyclonal Antibody-FITC conjugate
Storage	-20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>95%; Protein G purified
Buffer	10 mM Phosphate Buffered Saline, pH 7.4