

Product Datasheet

Murine IL-15 ELISA Kit 660.090

Product Name	Murine IL-15 ELISA Kit
Catalog Number	660.090
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	6.25-400pg/ml
Sensitivity	2.0pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant murine IL-15
Sample Types	Cell culture supernatant and serum tested, other body fluids suitable