

Product Datasheet

Alpha-fetoprotein (AFP) Antibody B01-94-02M-P

Product Name	Alpha-fetoprotein (AFP) Antibody
Catalog Number	B01-94-02M-P
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	AFP
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	0.05M Phosphate pH 7.15 with 0.85% Sodium Chloride
Specificity	AFP
Preservatives	0.1% Sodium Azide