

Product Datasheet

Methamphetamine Antibody C03-99-31M-P

Product Name	Methamphetamine Antibody
Catalog Number	C03-99-31M-P
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Methamphetamine
Matched Pairs	P81-99-31A (capture) or P91-99-31H (detection)
Preservatives	0.05% Sodium Azide