

# Product Datasheet

## Murine IL-17F ELISA Kit 660.100

Product Name	Murine IL-17F ELISA Kit
Catalog Number	660.100
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	62.5-4000pg/ml
Sensitivity	10pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant murine IL-17F
Sample Types	Cell culture supernatant and serum tested, other body fluids suitable