

Product Datasheet

Cytomegalovirus (CMV-gB) Antibody CA005

Product Name	Cytomegalovirus (CMV-gB) Antibody
Catalog Number	CA005
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	WB, ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2–8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	> 95 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	gB